SIDE ACCESS RACKS

Side Opening Horizontal units also known as side access racks, fit onto shelves of upright freezers to hold cryoboxes or microtiter plates. Side access racks feature a handle on one end to remove it from the freezer, these racks are best suited for longer term storage. Side access racks are manufactured from aluminium or stainless steel. They can be any size to meet the requirements of the user and have optional locking rods to hold items in place with a handle fitted to the top to make the unit easier to transport.

Labmode also offers a very popular adjustable side access rack which offers the user flexibility to store different size boxes. Using the unique design of this rack the user can set the removable shelves to store either $16 \times 2''$ boxes, $12 \times 3''$ boxes or $8 \times 4''$ boxes. Part Number SOR002

Dimensions (mm)	Box Size / Boxes Per Rack	Part Number
295H x 145W x 290D	4" (100mm) / 4	SOR455
295H x 145W x 290D	3" (75mm) / 6	SOR453
295H x 145W x 290D	2" (50mm) / 10	SOR110
420H x 140W x 560D	4" (100mm) / 12	SOR304
420H x 140W x 560D	3" (75mm) / 20	SOR368
420H x 140W x 560D	2" (50mm) / 28	SOR090
235H x 141W x 560D	2" (50mm) / 16	SOR196
235H x 141W x 560D	3" (75mm) / 12	SOR197
290H x 140W x 560D	2" (50mm) / 20	SOR032
340H x 140W x 560D	2" (50mm) / 24	SOR033
230H x 140W x 560D	2" (50mm) / 16	SOR196

Also available in stainless steel.

BESPOKE STORAGE

Should you require a bespoke front or side access rack please visit our website (www.labmode.co.uk) to fill in a quotation request. The information that will be required is the freezer make and model, maximum internal size of the shelf and the size of the item you require racking.



Side access rack to store 2" box SOR196



Side access rack to store microtiter plates SOR004



Side access rack to store various height boxes SOR002B

